

AED Maintenance Checklist

Hester / Funk /	′ 11 Apr 2025				Complete	
Score	0 / 0 (0%)	Flagged items	2	Actions	1	
First Name:					Hester	
Last Name:						
Role in the Company:					Safety Officer - APAC	
Is the AED in place?					Yes	
In this inspection, did you notice that the AED indicator or status light was flashing, or showing "rescue ready" or "operational"?					Yes	
Did you notice any cracks in the exterior of the device, or its sockets?					Yes	
Is the battery within the expiry date as shown on the battery or packaging?					Yes	
In this inspection, did you change the battery?					No	
Battery is okay. But also, we don't have extra battery in storage.						
To do Priority: High Due: 18.04.2025 15:13 PST Created by: SafetyCulture Staff						
Restock AED batteries						
Is there a set of sealed AED pads that are within the expiry date?					Yes	
In this inspection, did you change the pads?					Yes	
Check the following that apply:						
Prep pack (e.g., scissors, razor) present					v	
AED signage in place					v	
AED is visible					✓	
AED is accessible					✓	

Upload a short 5 second video of Status Indicator

aed-status-check-vid.mp4

Upload a photo of pad expiry date

Expiry date no longer readable on AED but still on box.



Photo 1

Upload a photo of battery expiry date

Expiry date no longer readable on battery but still on box. Came with AED pack as a set.



Photo 2

Do you want to make any notes on this maintenance log?

N/A

Signature:

Hester Funk 14.04.2025 10:41 PST

Date: 11.04.2025

Flagged items & Actions	2 flagged, 1 action					
Flagged items	2 flagged, 1 action					
AED Maintenance Log Did you notice any cracks in the exterior of the device, or its sockets?	Yes					
AED Maintenance Log In this inspection, did you change the battery?	No					
Battery is okay. But also, we don't have extra battery in storage.						
To do Priority: High Due: 18.04.2025 15:13 PST Created by: SafetyCulture Staff						
Restock AED batteries						
Other actions	0 actions					

Media summary





Photo 1 Photo 2

File summary

aed-status-check-vid.mp4